

Ordering Information

Endoscopic Approach					Percutaneous Approach										
Code	Stent		Delivery		Code	Stent		Delivery							
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)						
BD0604	6	4	8	180	TD0604	6	4	8	50						
BD0605		5			TD0605		5								
BD0606		6			TD0606		6								
BD0607		7			TD0607		7								
BD0608		8			TD0608		8								
BD0609		9			TD0609		9								
BD0610		10			TD0610		10								
BD0612		12			TD0612		12								
BD0804		8			4		8			180	TD0804	8	4	8	50
BD0805					5						TD0805		5		
BD0806	6		TD0806	6											
BD0807	7		TD0807	7											
BD0808	8		TD0808	8											
BD0809	9		TD0809	9											
BD0810	10		TD0810	10											
BD0812	12		TD0812	12											
BD1004	10	4	8	180	TD1004	10	4	8	50						
BD1005		5			TD1005		5								
BD1006		6			TD1006		6								
BD1007		7			TD1007		7								
BD1008		8			TD1008		8								
BD1009		9			TD1009		9								
BD1010		10			TD1010		10								
BD1012		12			TD1012		12								