



DeltaFlo

Orthopaedic suction set with filter



Distributed By

 **Instrumed Surgical**

Electro/Medical Sales and Service
A Div. of 707284 Ont. Inc.
2180 DUNWIN DRIVE, UNIT #5 & #6, MISSISSAUGA, ONTARIO, L5L 5M8
Customer Service #(800)-667-5653 ** Ph.#(905)-820-0902 ** Fax#(905)-820-0983
** Email: instrumed@instrumed surgical.com

More efficient More flexible More cost-effective



Removes debris and bone fragments without blocking

The large amounts of debris produced during orthopaedic surgery (bone fragments, chips, cement, marrow and blood clots) can often block conventional narrow-gauge suction sets. Now, thanks to DeltaFlo it's easy to keep the operative area visible and clean. The special filter built into the handle automatically traps solids and can be exchanged in a matter of seconds when filled. Not only does DeltaFlo minimize contact with body fluids, its effective non-block design helps shorten the procedure which in turn decreases the risk of damage to the tissue, infection and post-operative complications. Additionally, DeltaFlo significantly reduces damage to active red blood cells, which makes it ideal for collection of blood.

- Patented filter system prevents inconvenient and time-consuming blockages
- Increased filter capacity makes DeltaFlo ideal for hip arthroplasty and other procedures that produce a large amount of debris
- Ideal for peri-operative blood collection
- Also well-suited for treatment of aortic aneurysms, spinal surgery, and abdominal/urological procedures including removal of kidney stone fragments

Each DeltaFlo package includes an extra filter for quick exchange if the original filter becomes filled with debris.



Relieving eyes reduce the risk of trauma when working around soft tissue. The suction tip can be reversed for intra-femoral suction to remove fragments deep inside the femur.



- Non-block design decreases risk of damage to tissue and post-operative infection and complications
- Filter exchange takes only a few seconds
- Reversible suction tip provides extra length for intra-femoral suction and eliminates need for additional equipment
- Non-slip handle ensures total control, even with wet hands



The patented cross-flow filter allows fluids to pass through with minimal damage, yet traps debris that would clog conventional suction sets.

Material specifications

Handle:	Styrene Acrylonitrile
Suction tip:	Polyvinyl Chloride
Filter:	Styrene Acrylonitrile
Connecting tube:	Polyvinyl Chloride
Connector:	Polyvinyl Chloride
Double peel-packed Sterilized by ethylene oxide	

OP-Flex DeltaFlo

	Length of suction tip	Inner/outer diameter of suction tip	Length of tubing	CH/FG of connecting tube	REF
Suction set	230 mm	6 / 8 mm	2500 mm	24	D202
	150 mm	6 / 8 mm	2500 mm	24	D212
Suction handle	230 mm	6 / 8 mm	-	-	D102
	150 mm	6 / 8 mm	-	-	D112

All Unomedical DeltaFlo Products are Latex Free

Instrumented Product /Order Code for the complete set (handle, 2 filters, 3 tips) 33094182

For further information and availability please contact:

Unomedical Ltd.
Stroudwater Business Park,
Stonehouse, GL10 3SX.
Great Britain

Tel +44 (0)1453 827696
Fax +44 (0)1453 827861
E-mail bue@unomedical.com
www.unomedical.com

Stroud 0831U11.03

Distributed By



Electro/Medical Sales and Service

A Div. of 707284 Ont. Inc.

2180 DUNWIN DRIVE, UNIT #5 & #6, MISSISSAUGA, ONTARIO, L5L 5M8

Customer Service #(800)-667-5653 ** Ph.#(905)-820-0902 ** Fax#(905)-820-0983

** Email: instrumed@instrumedurgical.com