




Esophageal Stent



HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

Niti-S™ Esophageal Stent

for benign and malignant esophageal strictures

Features

- Fixed cell with braided construction
 - High flexibility and optimal radial force
 - Both head ends (8mm larger than trunk) help to minimize migration
- Silicone covering and soft round ends
 - Reduce tissue ingrowth and hyperplasia reaction
- Visible green suture for easy removal
- Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle

Ordering Information

| Distal Release | | | | | |
|---|---------|---------------|-------------------------|--------------|--------------------|
| Product | Code | Stent | | Delivery | |
| | | Diameter (mm) | Length (cm) | Profile (Fr) | Usable length (cm) |
| Fully covered | ES**_F | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | 70 |
| | | 20, 22 | | 20 | |
| | | 24, 28 | | 22 | |
| Both ends bare | ES**_B | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | |
| | | 20, 22 | | 20 | |
| | | 24, 28 | | 22 | |
| Fully covered / Removal string on both ends | ES**_F2 | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | |
| | | 20, 22 | | 20 | |
| | | 24, 28 | | 22 | |

* Available in USA (The sizes marked in blue)

| Proximal Release | | | | | |
|---|----------|---------------|-------------------------|--------------|--------------------|
| Product | Code | Stent | | Delivery | |
| | | Diameter (mm) | Length (cm) | Profile (Fr) | Usable length (cm) |
| Fully covered | ESP**_F | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | 70 |
| | | 20, 22 | | 20 | |
| | | 24, 28 | | 22 | |
| Both ends bare | ESP**_B | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | |
| | | 20, 22 | | 20 | |
| | | 24, 28 | | 22 | |
| Fully covered / Removal string on both ends | ESP**_F2 | 16, 18 | 6, 8, 10, 12, 14, 15 | 16 | |
| | | 20, 22 | | 20 | |
| | | 24 | | 22 | |

Released Article

* Fully covered, retrievable self-expanding metal stents (Niti-s) in palliation of malignant dysphagia : Long-term results of a prospective study by Sung Jun Choi et al [Scand J Gastroenterol. 2011 Jul;46(7-8):875-80]