




# CERVICAL™ Esophageal Stent



# Niti-S™ CERVICAL™ Esophageal Stent

for upper esophageal strictures

## Features

- Short proximal head design prevents damage of vocal cords in cases of stent placement close to the upper esophageal sphincter
- Silicone covering: Reduce the risk of tumor ingrowth
- Visible green suture for easy removal
- Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle

## Ordering Information

Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (mm)	Profile (Fr)	Usable Length (cm)
ES1606FV	16	6	18	70
ES1608FV		8		
ES1610FV		10		
ES1612FV		12		
ES1614FV		14		
ES1615FV		15		
ES1806FV	18	6	20	
ES1808FV		8		
ES1810FV		10		
ES1812FV		12		
ES1814FV		14		
ES1815FV		15		
ES2006FV	20	6	20	
ES2008FV		8		
ES2010FV		10		
ES2012FV		12		
ES2014FV		14		
ES2015FV		15		
ES2206FV	22	6	22	
ES2208FV		8		
ES2210FV		10		
ES2212FV		12		
ES2214FV		14		
ES2215FV		15		

Proximal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (mm)	Profile (Fr)	Usable Length (cm)
ESP1606FV	16	6	18	70
ESP1608FV		8		
ESP1610FV		10		
ESP1612FV		12		
ESP1614FV		14		
ESP1615FV		15		
ESP1806FV	18	6	20	
ESP1808FV		8		
ESP1810FV		10		
ESP1812FV		12		
ESP1814FV		14		
ESP1815FV		15		
ESP2006FV	20	6	20	
ESP2008FV		8		
ESP2010FV		10		
ESP2012FV		12		
ESP2014FV		14		
ESP2015FV		15		
ESP2206FV	22	6	22	
ESP2208FV		8		
ESP2210FV		10		
ESP2212FV		12		
ESP2214FV		14		
ESP2215FV		15		

Distal Release				
► Removal string on both ends				
Code	Stent		Delivery	
	Diameter (mm)	Length (mm)	Profile (Fr)	Usable Length (cm)
ES1606FVR2	16	6	16	70
ES1608FVR2		8		
ES1610FVR2		10		
ES1612FVR2		12		
ES1614FVR2		14		
ES1615FVR2		15		
ES1806FVR2	18	6	20	
ES1808FVR2		8		
ES1810FVR2		10		
ES1812FVR2		12		
ES1814FVR2		14		
ES1815FVR2		15		
ES2006FVR2	20	6	20	
ES2008FVR2		8		
ES2010FVR2		10		
ES2012FVR2		12		
ES2014FVR2		14		
ES2015FVR2		15		
ES2206FVR2	22	6	22	
ES2208FVR2		8		
ES2210FVR2		10		
ES2212FVR2		12		
ES2214FVR2		14		
ES2215FVR2		15		
ES2406FVR2	24	6	22	
ES2408FVR2		8		
ES2410FVR2		10		
ES2412FVR2		12		
ES2414FVR2		14		
ES2415FVR2		15		

Proximal Release				
► Removal string on both ends				
Code	Stent		Delivery	
	Diameter (mm)	Length (mm)	Profile (Fr)	Usable Length (cm)
ESP1606FVR2	16	6	16	70
ESP1608FVR2		8		
ESP1610FVR2		10		
ESP1612FVR2		12		
ESP1614FVR2		14		
ESP1615FVR2		15		
ESP1806FVR2	18	6	20	
ESP1808FVR2		8		
ESP1810FVR2		10		
ESP1812FVR2		12		
ESP1814FVR2		14		
ESP1815FVR2		15		
ESP2006FVR2	20	6	20	
ESP2008FVR2		8		
ESP2010FVR2		10		
ESP2012FVR2		12		
ESP2014FVR2		14		
ESP2015FVR2		15		
ESP2206FVR2	22	6	22	
ESP2208FVR2		8		
ESP2210FVR2		10		
ESP2212FVR2		12		
ESP2214FVR2		14		
ESP2215FVR2		15		
ESP2406FVR2	24	6	22	
ESP2408FVR2		8		
ESP2410FVR2		10		
ESP2412FVR2		12		
ESP2414FVR2		14		
ESP2415FVR2		15		