

# H. pylori Tests

## Presumptive Identification of H. pylori



**HpFast™**  
Results in 24 hours

- Both HpFast and HpOne provide a positive control kit, a breakaway biopsy pick and spatula
- Easy-to-read positive indicator gives you a quick, timely result
- No need for unsafe use of a needle or syringe

**HpOne™**  
Results in 1 hour



HpOne™ and HpFast™ kits include everything needed to conduct a quick and accurate urease test for the presumptive identification of Helicobacter. The easy-to-read indicator make Hp tests by GI Supply the clear choice in H. pylori testing.



## HpFast™

The agar gel formula of this test provides accurate H. pylori status results in **24 HOURS** or less.

## HpOne™

This test provides accurate, cost-effective H. pylori status results in **ONE HOUR** or less. Simple. HpOne liquid reagent included.

## Ordering Information

Item No.	Title	Quantity
GIS-22M	HpFast	Box of 50
GIS-22D	HpFast	Box of 20
GIS-64M	HpOne	Box of 50
GIS-64D	HpOne	Box of 20
GHP-U	Urease	1 ML bottle

HpFast comes with positive control (Urease) and breakaway biopsy pick.

HpOne comes with positive control (Urease), breakaway biopsy pick and Reagent.

- CLIA waived (Clinical Laboratory Improvement Amendments)
- Reimbursed under CPT Code 87077QW -  
Presumptive identification of Helicobacter pylori in gastric mucosal biopsies



An ISO 13485 Certified Company